



OR 4Vet

Multimodal image diagnostics  
that inspire.

# Leonardo DR mini III

## The premium X-ray case

with large diagnostic monitor and extremely long battery runtime



 OR Technology

The expert by your side.



# The gold standard for equine and mixed practices



The case can be comfortably transported by the **handle** or with the **carrying strap**

Professional software for detailed and complete **documentation of X-ray services for horse purchase examinations** (for Germany only) according to the „X-ray Guide 2018“ (optional)

**Intelligent power saving mode** - fold down the monitor for standby and continue immediately within a certain time frame. System shuts down automatically if not re-opened.



The integrated monitor is brought into the working position with a flick of the wrist

Robust, **extremely durable suitcase** made of the latest high-tech composite material with shock-absorbing edge reinforcement, splash-proof (IPX4)



Optional **wireless remote control** of the system via your smartphone using an app (see video link on the back)

**Unbeatable image quality** through automatic, intelligent image processing

**Very low case weight** of only 9.5 kg with large Full HD touchscreen monitor (without detector)



Sufficient space for a **25 x 30 cm** or **35 x 43 cm** detector incl. protection box



Easy removal of the X-ray detector from the **padded transport compartment**

**Anti-glare 21.5" (54.6 cm) Full HD touchscreen monitor** for large image display and viewer comfort diagnosis



**Intuitive operation** via *dicomPACS® DX-R* acquisition software with **integrated X-ray positioning guide** - separated according to animal species (horse, dog, cat, etc.)



Power can also be supplied with an **external power pack** without the need for batteries\*

\* suitable for emergency mode only

**Status display** of the overall system in 4 colours for direct visualisation of the workflow

Standard **PC keyboard**

**USB connection**



**Replacing the batteries during operation**, no need to shut down the system

Optional **accessory bag** for safe transport of spare batteries, chargers etc.



**Special vet measuring** and macros available for recurring examination procedures, e.g. open angle for laminitis, purchase examination etc.

System enables continuous, **cordless work for 8 hours with up to 500 X-ray exposures** (in double-battery operation; 18 Volt / 5 Ah)



## Leonardo DR mini III

# Almost indestructible 25 x 33 cm CsI detector without glass carrier

**World's first CsI X-ray detector with flexible substrate - specially developed for imaging in horses**

The compact, robust X-ray detector XenOR 33-100W-flex with flexible carrier material (no glass) incl. CsI scintillator was specially developed for use in equine medicine. With a pixel density of 100  $\mu\text{m}$ , the X-ray detector guarantees the highest resolution for razor sharp and detailed images. It weighs only approx. 2 kg, withstands drop heights of 120 cm without damage and has an almost waterproof housing.

The narrow edges of approx. 7 mm (on two sides) ensure optimum hoof, knee and tooth images. The battery charge lasts up to 8.5 hours and the cycle time between images is < 3.5 sec – this guarantees fast work.

- Automatic synchronisation of detector and generator via AED - no need to intervene in the X-ray system or adjust the system or cable connections
- Battery charge lasts for approx. 8.5 hours
- Narrow edges of only approx. 7 mm on 2 sides
- Degree of protection IP67
- Active area:  
252.8 mm x 331.8 mm [10" x 13"]
- Pixel density 100  $\mu\text{m}$
- Pixel matrix 3318 x 2528
- Weight approx. 2 kg
- Dimensions (W x H x D): 269 x 362 x 14.5 mm





## Leonardo DR mini III

# Long running times through intelligent battery concept

### Up to 500 shots without charging

A completely new concept for the power supply has been developed for the Leonardo DR mini III, allowing you to choose the batteries. Depending on intended use, you can X-ray for 3.5, 8.5 or unlimited hours.

Batteries can be changed during operation.

→ **Advantage:** No shutdown of the system is necessary.

Batteries are not discharged when the case is switched off.

→ **Advantage:** Case is ready for immediate use, even if it has not been used for a long time.

Unit automatically switches to stand-by mode when the lid is closed.

→ **Advantage:** After changing the location, the computer does not have to be restarted.

The standard rechargeable batteries can be bought in almost any DIY store worldwide.

→ **Advantage:**

- no expensive special parts with complex logistics
- choice of large or small batteries enable user to save weight or extend running time as required.
- With low financial investment, several batterysets can be used and, depending on requirements, an unlimited runtime can be achieved

Unlimited  
power  
possible



**Battery set (small)**  
18 Volt/2 Ah for a runtime of approx. 3.5 hours and approx. 200 X-ray exposures (0.35 kg/battery)



**Battery set (large)**  
18 Volt/5 Ah for a running time of approx. 8.5 hours and approx. 500 X-ray exposures (0.62 kg/battery)

# Advantages of the professional acquisition software

Modern graphical user interface, **touchscreen** operation

Capture of patient data via **DICOM Worklist**, **BDT/GDT**, **HL7** or other protocols – data may also be captured manually

Use of **DICOM Procedure Codes** for the transfer of all relevant examination data directly from the connected patient management system (HIS/RIS)

**Freely configurable** body parts with more than **200 projections** and numerous possible adjustments

Safe and fast **registration of emergency patients** allowing the user to **switch between examinations** of a patient, for instance to avoid having to re-position the patient frequently

Allows the user to **subsequently add images** to an examination, even after that examination has already been completed

Special tools for veterinary medicine, such as an extra dialogue box for patient and owner data, **automatic hip dysplasia measuring (AI-based)**, **MMP**, **TTA**, **TPLO** and **vertebral heart score (Buchanan)**, **special image filters** and much more

Fully **integrated radiographic positioning guide** for each examination incl. comprehensive notes, photos and correct X-ray images

Suited for use in dental medicine, e.g. with a digital dental sensor or a dental CR system





## Leonardo DR mini III

# Standard components of the X-ray case

### Compact, lightweight X-ray suitcase

- Extremely durable case made of the latest high-tech composite material with shock-absorbing edge reinforcement, splash-proof (IPX4), lockable
- Total weight: approx. 9.5 kg (includes complete X-ray case, monitor and and electronics, plus X-ray detector)
- Anti-glare 21.5" full HD touchscreen monitor and standard PC keyboard
- High-performance PC with current Windows version
- Carrying strap for comfortable transport
- LED status display of the entire system in 4 colours
- Padded transport compartment for one detector, max. 14" x 17"
- Dimensions: approx. 545 x 515 x 194 mm (W x D x H)
- Connections: 1 x USB, 1 x LAN, 1 x WLAN, 1 x Bluetooth, 1 x DC socket



### CsI X-ray detector flex 25 x 33 cm (13" x 10") wireless

Wireless X-ray!

World's first CsI X-ray detector with flexible substrate - almost indestructible - specially developed for imaging horses

or

### CsI X-ray detector flex 35 x 43 cm (14" x 17") wireless

Wireless X-ray!

World's first CsI X-ray detector with flexible substrate - almost indestructible - specially developed for imaging horses

Various types of X-ray detectors on request



### Leonardo DR mini III software package

with dicomPACS®DX-R, professional console software with modern graphical user interface including basic software package and integrated radiographic positioning guide for each examination, inclusive:

- dicomPACS®DX-R DICOM Send SCU
- dicomPACS®DX-R DICOM Patient CD
- dicomPACS®DX-R Cognition Optimised Processing



### Battery set (optionally small or large)

- Battery set (small) 18 Volt/2 Ah for a runtime of approx. 3.5 hours and approx. 200 X-ray exposures (0.35 kg/battery)
- Battery set (large) 18 Volt/5 Ah for a running time of approx. 8.5 hours and approx. 500 X-ray exposures (0.62 kg/battery)



## Optional components for the Leonardo DR mini III system

### Wireless remote control

Possibility of wireless control of the digital X-ray system with the app "dicomPACS®DX-R remote control" via smartphone or tablet incl. display of the worklist, preview of the captured image to check the shot etc. - directly from the app a shot can be repeated or a new one can be planned



Video for wireless remote control of the system via your smartphone (via app) under [www.or-technology.com/DX-R-remote-control-video-vet-en](http://www.or-technology.com/DX-R-remote-control-video-vet-en)



### Protection case 1417

Dimensions 47.9 x 46.7 x 2.54 cm (W x H x D), weight approx. 1.35 kg

### Accessories bag

for transport and storage of spare batteries, chargers, cables, power supply etc.

### Amadeo P high-frequency X-ray units

Portable monobloc X-ray units for high-quality X-rays: Low weight and simple operation guarantee versatile use in small animal practices and equine clinics







# OR 4Vet

Multimodal image diagnostics  
that inspire.

**Everything from a single source:** PACS | X-ray | Ultrasound | CT  
OR 4Vet. Because image diagnostics only needs one partner.

Headquarters:

OR Technology (Oehm und Rehbein GmbH), Neptunallee 7c, 18057 Rostock, Germany

Tel. +49 381 36 600 500, Fax +49 381 36 600 555

[www.or-technology.com](http://www.or-technology.com), [info@or-technology.com](mailto:info@or-technology.com)

OR Technology UK: Celtic SMR Ltd., Frederick House, Hayston View, Johnston  
Haverfordwest, Pembrokeshire SA62 3AQ, United Kingdom

[www.celticsmr.co.uk](http://www.celticsmr.co.uk), [sales@celticsmr.co.uk](mailto:sales@celticsmr.co.uk)

Info hotline: +49 381 36 600 600

[www.or4vet.com](http://www.or4vet.com)



More information about  
Leonardo DR mini III systems  
[www.oehm-rehbein.de](http://www.oehm-rehbein.de)

Specifications subject to revision without notice

The editor strives to impart correct and up to date information. The provided specifications are based on current knowledge and are subject to revision without notice. This brochure is subject to correction. The editor assumes no responsibility for the information being up to date, correct and complete.

All furnished logos, pictures and graphics are property of the particular company and subject to copyright of the licensor. Use, dissemination, distribution or copying of the pictures, logos or text compiled or processed by the editor is subject to our written consent. All rights reserved.

