


Benefits

Outstandingly sophisticated image diagnosis

- Completely integrated **dicomPACS® Viewer for image diagnosis**, further processing and storage of images in an SQL database incl. image manipulations, export options, layout adjustments, freely configurable user interface and much more



- Stepless **zoom, PAN, magnifyer, ROI, crop, rotate, mirror** etc.
- Insertion of **image annotations**, e.g. free texts, arrows, ellipses etc.
- Measuring** of distances, angles, areas and density
- Adjustment of window/level options and **gamma correction**, sharpening filters, noise suppression
- Many additional functions such as calculation of **Cobb's angle, HD measurements, pelvic obliquity measurements, integrated capturing of diagnostic reports** etc.
- Printing of images both on Windows printers and laser imagers via **DICOM Basic Print**
- Creation of **DICOM patient CDs** with free **WEB viewer**
- Export of images** to JPEG, TIFF, BMP and DICOM formats
- Easily upgradable to the professional, **integrated image management system (PACS)**