



Product·Information

Mobile X-ray table with integrated holder for cassette and grid

Brief description

- Integrated holder for cassette and grid under the table top along the table
- Ideal for use in combination with our wall mounts
- Table frame: anodised aluminium
- Table top: radiolucent white synthetic
- Easy to clean and disinfectant proof
- Optional: adjustable positioning aide for dysplasia coxae
- Optional: wall bracket and/or sprung wall bracket for X-ray unit
- Dimensions:
 Height: 78.0 cm
 Length: 155.5 cm
 Depth: 67.5 cm

- **Cassette and grid holder under the table top, slides along the table**



Mobile X-ray table with integrated holder for cassette and grid



Specifications subject to revision without notice

The editor strives to impart correct and up to date information. The provided specifications are based on current knowledge and are subject to revision without notice. This brochure is subject to correction. The editor assumes no responsibility for the information being up to date, correct and complete.

All furnished logos, pictures and graphics are property of the particular company and subject to copyright of the licensor. Use, dissemination, distribution or copying of the pictures, logos or text compiled or processed by the editor is subject to our written consent. All rights reserved.

Stemp of distribution partner