

# Product·Information

## **dicomPACS® DX-R 19" touch screen workstation** [Mini PC with 19" touch screen monitor]

### Scope of delivery

- Mini PC
- 19" touch screen monitor
- Optical mouse and keyboard
- MS Windows Professional
- Hardware stress test
- **dicomPACS® DX-R** software comes pre-installed

### Specification

#### Mini PC

- Intel® Core processor
- 4 GB RAM DDR
- 2 x 500 GB hard disks configured as mirror (RAID1)
- DVD+/-RW drive
- 2 x 1000 Mbit network card

#### 19" touch screen monitor color

ProLite T1931SR-1 is based on the latest 5 Wire Resistive technology. Price / performance ratio is ideal for the budget conscious user without compromising on quality and specifications.

- 19" diagonal (Pixel pitch h x w 0.294 x 0.294 mm)
- Resistive touch technology
- 80% light transmittance
- Display area h x w 301 x 376 mm; 11.9" x 14.8"
- 5 ms response time
- Contrast 900 : 1 with touch panel
- 200 cd/m<sup>2</sup> brightness with touch panel
- Viewing zone horizontal/vertical: 170°/160°;  
right/left: 85°/85°; up/down: 80°/80°
- 16.7 million displayable colours



# dicomPACS® DX-R 19" touch screen workstation

## Dimensions:

### Monitor



### Desktop



### Mini-PC

*[Please plan for sufficient space for plug connectors and ventilation!]*



### Specifications subject to revision without notice

The editor strives to impart correct and up to date information. The provided specifications are based on current knowledge and are subject to revision without notice. This brochure is subject to correction. The editor assumes no responsibility for the information being up to date, correct and complete.

All furnished logos, pictures and graphics are property of the particular company and subject to copyright of the licensor. Use, dissemination, distribution or copying of the pictures, logos or text compiled or processed by the editor is subject to our written consent. All rights reserved.

Distribution partner