

References

Veterinary Hospital "Am Bökelberg"
Dr. med. vet. Arnd Stelljes,
41063 Mönchengladbach, Germany



Dr Arnd Stelljes: "Our veterinary hospital has set itself the task of being able to offer - by virtue of specialisation - a comprehensive, state of the art small animal veterinary service to its clients. All our diagnostic practices are continually adapted to the latest scientific developments.

Our clinic has two areas of particular specialisation: orthopaedics and surgery. To minimise X-ray exposure and ensure the quick and reliable management of examination results, we have installed an advanced digital diagnostic system.

An AGFA CR 25.0 with a special veterinarian user interface by GIERTH X-Ray international was installed, while an Esaote Vet-MR system as well as an ultrasound medicine from GE was also integrated into our practice network via the PACS software.



The image processing software **dicomPACS®vet** by OR Technology could be linked directly and easily to our patient management system (Vetera by GP Software). Now all X-rays, MRI and ultrasound images can be easily displayed and processed at our six workstations in the treatment rooms.

The direct interface of the image processing software with our patient management system offers us preview images from all linked modalities in our Vetera program. Images can now also be directly called up from the card index. After initial network problems affected loading times, the image transfer rate is now also highly satisfactory.

In keeping with the saying "To each his own" we opted for the integration of a pure PACS software with a proven patient management system.

In my view **dicomPACS®vet** is a PACS programme that will fully meet our needs, now and in the future".



Innovative image management solutions

Nowadays, diagnostic investigation with the help of digital image solutions plays an important role in the daily work of each modern practice or hospital. As the use of these modern technologies will continue to increase during the next few years, we would like to give you an understanding of our PACS software **dicomPACS®vet**.

dicomPACS®vet is an innovative and intelligent high tech solution for image processing as well as practice and hospital management; with the help of **dicomPACS®vet** a paperless daily working routine can become reality.

All types of images (X-rays, CT, MRI, digital camera, ultrasound) as well as all kinds of documents (doctors' letters, diagnoses, recovery processes, faxes) can be acquired, diagnosed and then filed in the digital patient file with **dicomPACS®vet**, only one mouse click away.

We guarantee quick access to all data under the highest international security standards with our well designed archive and backup solutions. Moreover, **dicomPACS®vet** can be integrated effortlessly with management systems via HL7 or BDT/GDT communication or using an individual solution.



OR Technology (Oehm und Rehbein GmbH)
18057 Rostock, Germany, Waldemarstr. 20 g/h
Info-Hotline: +49 (0)381 - 20 36 126
www.or-technology.com, info@or-technology.com

[Stamp of distribution partner]

Further information is available under www.or-technology.com