

Compact, mobile X-ray system **Leonardo DR Go** e.g. for tuberculosis control

AI-based
thoracic screening
for automatic
tuberculosis
detection
– optional –

- Easy and secure transport of the complete digital X-ray system in the indestructible Peli™ Storm Case™, including:

- Generator
- X-ray detector
- Laptop including Leonardo software package
- Generator stand
- Detector mounts
- All accessories such as charger and charging cable

- High quality and robustness, very durable and suited for harsh environmental conditions

- X-ray equipment is securely integrated into customised inserts made of high-quality polyurethane foam in the case

- Standard locking system

- Polyurethane wheels that can withstand rough surfaces

- 2-stage extendable handle

- Dustproof and splash-proof
(certified according to IP 67)

Compact, transportable X-ray system Leonardo DR Go

Parameter

General	
External dimensions (LxWxH)	63.2 × 60.2 × 33.3 cm
Colour	black
Weight including accessories	ca. 28 kg

Optional	
AI-enabled notebook with integrated tool for AI-based chest screening for automatic tuberculosis detection	

Suitcase features	
6 push-pull fasteners, telescopic handles	
2 double-layered, soft-touch handles	
2 padlock-lockable snap locks	
Sturdy hinges, O-ring seal	

Integrated accessories	
X-ray generator: Amadeo P-90/20VB – lightweight, up to 300 images	
X-ray detector: XenOR 35-100WF for use up to 8 hours	
Acquisition console: <i>dicomPACS® DX-R</i> Notebook 16"	
Detector mount: VersariX – flexible and lightweight (400 g)	
Generator stand: height-adjustable carbon tripod	
Generator and detector accessories	



Please visit www.or-technology.com
for more details

Specifications subject to revision without notice

The editor strives to impart correct and up to date information. The provided specifications are based on current knowledge and are subject to revision without notice. This brochure is subject to correction. The editor assumes no responsibility for the information being up to date, correct and complete. All furnished logos, pictures and graphics are property of the particular company and subject to copyright of the licensor. Use, dissemination, distribution or copying of the pictures, logos or text compiled or processed by the editor is subject to our written consent. All rights reserved.

Headquarters: OR Technology (Oehm und Rehbein GmbH), 18057 Rostock, Germany, Neptunallee 7c
Info-Hotline: +49 381 36 600 600, info@or-technology.com, www.or-technology.com
OR Technology UK: Celtic SMR Ltd., Frederick House, Hayston View, Johnston, Haverfordwest, Pembrokeshire SA62 3AQ, United Kingdom, www.celticsmr.co.uk, sales@celticsmr.co.uk

