

Reference

Direct digital X-ray with **Amadeo** in the joint veterinary practice of Dr. Beate Fischer-Wackes and Dr. Jörg Wackes in Gundelfingen, Germany



Veterinary group practice

Dr. Jörg Wackes

Dr. Beate Fischer-Wackes

79194 Gundelfingen, Germany

Tel. +49 761-5950015

info@wackes-tieraerzte.de

www.wackes-tieraerzte.de

 Amadeo Systems

Reference

Dr. Beate Fischer-Wackes and Dr. Jörg Wackes offer all the services of a practice for small animals, including X-ray imaging, ultrasound and endoscopy. The focus is on the fields of dermatology, gold implantation, (endoscopic) surgery, dentistry and physiotherapy. Several hundred X-ray images are taken each quarter.

Since early 2011, the practice has been working with a direct digital X-ray system. The Amadeo DR system includes a patient table with a floating tabletop, the X-ray generator which is integrated in the X-ray table to save space, and a DR Tech Flaatz 560 detector.

“Before we purchased the new X-ray system we made enquiries amongst our colleagues, a number of whom recommended solutions by OR Technology”, says Dr. Jörg Wackes. “Our previous X-ray system was conventional. Using direct digital X-ray imaging with a detector made our workflow significantly faster”, he explains the advantages of the new X-ray system. “The image quality in full image mode is very good. Image processing and image weighting are available and useful, allowing us to obtain more information from one image.”

The *dicomPACS®DX-R* acquisition software is operated via a touch screen. The generated images are then evaluated and archived using the *dicomPACS®vet* image management system.

Further information is available under
www.or-technology.com



Dr. Jörg Wackes,
Specialist veterinarian for
small animals and pets



“Diagnostic evaluation is simple and up-to-date and the large number of *dicomPACS®vet* functions, such as annotations, grid displays, enlargement, zoom, light/dark contrast changes and image export are most helpful. Thanks to the uncomplicated archiving option we have fast access to images even at a later stage. The integration with our practice software Vetinf is working smoothly”, says Dr. Wackes.

Service and consultation are also important issues for veterinarians. “We are very happy with the technical support offered by OR Technology”, comments Dr. Wackes. “Any reported faults are processed by friendly telephone staff. In urgent situations, the problem is forwarded immediately or dealt with promptly. All in all, we can count on clear directions and solution-orientated help from OR Technology.”





The Amadeo X-ray system of the veterinary group practice Dres. Wackes was installed in 2011 and is no longer available. Below you will find our current model.

Space-saving
full system with rotating
monitor arm for use with
standard wall sockets



Amadeo

V *mini*

The Amadeo V-DR mini X-ray system has been developed especially for the requirements of veterinary practices and is particularly suited for small animal practices. The patient table is equipped with a floating tabletop and silent brakes. It allows for generous movement and can be operated from the front panel. The height-adjustable monitor arm can be rotated by 180°, enabling comfortable operation of the *dicomPACS®DX-R* acquisition and diagnostic software. The 230 V high-frequency generator is incorporated directly into the X-ray table, taking up minimal space. The compact build of the Amadeo V-DR mini makes it ideal for installation in confined spaces.

 OR Technology

www.or-technology.com

OR Technology (Oehm und Rehbein GmbH), 18057 Rostock, Germany, Neptunallee 7c
Tel. +49 381 36 600 500, Fax +49 381 36 600 555
www.or-technology.com, info@or-technology.com

Info-Hotline: +49 381 36 600 600



Further information
about the
Amadeo V mini
system is available here: