

Order No. AM1091

Plastic transport box for **Amadeo M mini** X-ray system

The PVC transport box offers optimal protection for transporting the **Amadeo M mini**. Loading the system into the box is easy thanks to the integrated ramp; there is no need to lift the system at all. The X-ray system can be transported in one piece, final assembly is not required.

- Well suited for quick relocation of the X-ray system
- Ribbed outer shell for extra sturdiness
- Developed for outdoor use: water repellent, shock proof, dust-tight
- High quality workmanship
- Foam padding inside the box encloses the X-ray system for a cushioning effect
- Highest transport safety due to robust design and special snap locks
- Loading via convenient built-in ramp
- Two trestle rolls and conveniently placed handholds make moving the box easy



Plastic transport box for **Amadeo M mini** X-ray system

Parameters


General	
Exterior measurements (L x W x H)	1560 x 870 x 810
Colour	white
Material	polyethylene
Empty weight	62.5 kg
Complete weight incl. X-ray system	141.5 kg (accessories not included)

Features


Transport box with built-in ramp
2 eyelets with elastic band to hold back the system while uprighting the system
8 front closures
14 handles
2 80mm trestle rolls

Further information

More information
about **Amadeo M mini**
you can be found here:
[www.or-technology.com/
AmadeoM-mini](http://www.or-technology.com/AmadeoM-mini)



Two sturdy trestle rolls for easy movement – even when loaded



Please visit www.or-technology.com
for more details

Specifications subject to revision without notice

The editor strives to impart correct and up to date information. The provided specifications are based on current knowledge and are subject to revision without notice. This brochure is subject to correction. The editor assumes no responsibility for the information being up to date, correct and complete. All furnished logos, pictures and graphics are property of the particular company and subject to copyright of the licensor. Use, dissemination, distribution or copying of the pictures, logos or text compiled or processed by the editor is subject to our written consent. All rights reserved.