

Order No. AM1090

# Multiplex wood transport box for the **Amadeo M mini** X-ray system

The multiplex transport box is very weatherproof and durable and offers optimal protection during the transport of the **Amadeo M mini**.

The box can be loaded without having to lift the X-ray system: The front can be folded down as a ramp and the system can be wheeled into the box.

The X-ray system can be transported in one piece, final assembly is not required.

- Multiplex plates with waterproof glue ensure stability and shape retention
- High quality workmanship
- Foam padding inside the box encloses the X-ray system for a cushioning effect
- Highest transport safety due to robust design and special snap locks
- Loading via convenient built-in ramp



# Multiplex wood transport box for the Amadeo M mini X-ray system

## Parameters

### General


Interior measurements (LxWxH)	990 x 8050 x 1700 mm
Material	wood (Multiplex)
Empty weight	170 kg
Complete weight incl. X-ray system	250 kg (accessories not included)

### Features

- Mounting systems with integrated foam padding
- 2 side restraint systems (have to be removed before loading and then fixed with two handy star knob screws)
- 6 safety latches in front, incl. lockable safety bolt (for padlock / combination lock)
- 2 side handholds
- 4 side block feet

### Further information

More information about **Amadeo M mini** you can be found here: [www.or-technology.com/AmadeoM-mini](http://www.or-technology.com/AmadeoM-mini)



Please visit [www.or-technology.com](http://www.or-technology.com) for more details

### Specifications subject to revision without notice

The editor strives to impart correct and up to date information. The provided specifications are based on current knowledge and are subject to revision without notice. This brochure is subject to correction. The editor assumes no responsibility for the information being up to date, correct and complete. All furnished logos, pictures and graphics are property of the particular company and subject to copyright of the licensor. Use, dissemination, distribution or copying of the pictures, logos or text compiled or processed by the editor is subject to our written consent. All rights reserved.