


Order No. AM 2200

Foldable X-ray table M-DR for mobile X-ray systems and X-ray machines

- Assembly and disassembly without tools in 2 minutes
- Weight only 48.5 kg*
- Transport case with castor wheels



**without attachments*

Foldable X-ray table M-DR

for mobile X-ray systems and X-ray machines



Advantages of the foldable X-ray table:

- Assembly and disassembly without tools:
2 persons, 2 minutes, 6 knurled screws
- Weight without attachments only 48.5 kg
- ESD castor wheels with breaks
- Movable infusion holder
- Movable mounting for lateral exposures
- Free positioning of the bucky tray under the entire table surface (detector can be rotated 90° for pelvic exposures)
- Detector fits into the bucky tray incl. protection case
Advantage: The detector is always protected and does not have to be removed from the case before each X-ray procedure
- Mounting for a removable grid
- Standard transport case with castor wheels incl. compartment for storage of detector and accessories
- Easy-to-clean table top



Grid optional



Please visit www.or-technology.com
for more details

Specifications subject to revision without notice

The editor strives to impart correct and up to date information. The provided specifications are based on current knowledge and are subject to revision without notice. This brochure is subject to correction. The editor assumes no responsibility for the information being up to date, correct and complete. All furnished logos, pictures and graphics are property of the particular company and subject to copyright of the licensor. Use, dissemination, distribution or copying of the pictures, logos or text compiled or processed by the editor is subject to our written consent. All rights reserved.

Headquarters: OR Technology (Oehm und Rehbein GmbH), 18057 Rostock, Germany, Neptunallee 7c

Info-Hotline: +49 381 36 600 600, info@or-technology.com, www.or-technology.com

OR Technology UK: Celtic SMR Ltd., Frederick House, Hayston View, Johnston, Haverfordwest, Pembrokeshire SA62 3AQ, United Kingdom, www.celticsmr.co.uk, sales@celticsmr.co.uk