

Wireless flat panel detector Varex-PerkinElmer XRPad2 4336

14 x 17 Zoll [35 x 43 cm]

XRpad2 4336 is a lightweight wireless flat panel detector with advanced features enabling digital radiography never seen before.

Featuring best-in-class 100 μm pixel size, direct deposition CsI scintillator or Gadox and excellent DQE, XRpad2 4336 enables high resolution imaging with reduced X-ray exposure. The design is lightweight, robust and ergonomic. Automatic Exposure Detection, on-board corrections, and wireless access point mode make system integration quick and simple.

- Dynamic mode with continuous output of up to 8 fps
- Onboard preview: The flat panel detector sends a preview immediately after image acquisition
- Acknowledge function: Buffering of acquired images
- High resolution 100 μm pixel pitch (5.0 lp/mm)
- Pixel array 3524 \times 4288
- Dimensions (W x H x D)
384 mm \times 460 mm \times 15.5 mm
- Weight approx. 3.2 kg (7.0 lbs)



Wireless flat panel detector

Varex-PerkinElmer XRPad2 4336

14 x 17 Zoll [35 x 43 cm]

Parameter

Sensor	
Scintillator	Direct deposition CsI:Tl or Gadox
Active Area	350 mm × 426 mm
Pixel Matrix	3524 × 4288
Pixel Pitch	100 μm

Electronics	
Amplifiers	Low noise ASICs with user selectable gains
ADC	16-bit
Image Transfer Time	Wired: 500 ms; Wireless: 3000 ms
On-board Memory	1 GB DDR3, 8 GB SDHC card

Mechanical	
Size	ISO 4090 for 35 cm × 43 cm (14" × 17") cassette size
Active Area	350 mm × 426 mm
External Dimensions (w/l/h)	384 mm × 460 mm × 15.5 mm
Weight	3.2 kg (7.0 lbs)
Housing	Aluminum frame with carbon-fiber entrance window

Communications	
Status Display	OLED display with Wi-Fi, LAN, battery, and sensor indicators
Wireless Data I/F	802.11n Wi-Fi standard at 5 GHz
Wired Data I/F	GigE, trigger and power via docking connector
X-ray I/F	Integrated X-ray trigger control Automatic Exposure Detection (AED)

Imaging Performance	
Typical DQE	75% (0 cy/mm), 60% (1 cy/mm), 40% (3 cy/mm) for RQA5
Typical MTF	70% (1 cy/mm), 40% (2 cy/mm), 15% (4 cy/mm) for RQA5
Limiting Resolution	5 cy/mm



Environmental	
Temperature	10 °C to 35 °C operating
Humidity	10% to 80% operating
Ingress Protection	IPX4 rated (splash water protection)

Accessories	
Battery	Rechargeable battery, 11.1 V
Battery Charger	External two bay charger 100 - 240 V, AC 50/60 Hz
Interface and Power Unit	Optional IPU-2 external power supply 100 - 240 V AC, GigE and X-ray I/F

Please visit www.or-technology.com for more details

Specifications subject to revision without notice

The editor strives to impart correct and up to date information. The provided specifications are based on current knowledge and are subject to revision without notice. This brochure is subject to correction. The editor assumes no responsibility for the information being up to date, correct and complete. All furnished logos, pictures and graphics are property of the particular company and subject to copyright of the licensor. Use, dissemination, distribution or copying of the pictures, logos or text compiled or processed by the editor is subject to our written consent. All rights reserved.



Battery charger with batteries