

# Reference

Upper Cervical (NUCCA) and Standard Chiropractic Tools,  
Dr. Jeffrey N. Scholten, „The Vital Posture™ Clinic“  
in Calgary, Canada



## Dr. Jeffrey N. Scholten

(President of Upper Cervical Health Centers Inc., Vice-President of both the National Upper Cervical Chiropractic Association (NUCCA) and the Upper Cervical Research Foundation (Canada), member of the board of directors of the International Chiropractic Association's Council on Upper Cervical Care, the Upper Cervical Research Foundation (UCRF)



# Reference



Our software for NUCCA (National Upper Cervical Chiropractic Association) and „upper cervical“-doctors offers all necessary tools to carry out NUCCA standards of practice and simplifies the workflow.

The system assists the users in all stages und gives instructions optionally. As a result, it is ensured that measurements are executed correctly at the right time. The software was developed and tested in closed cooperation by NUCCA-doctors from USA and Canada.

Doctor Jeff Scholten is a great advocate of the NUCCA software development and works with **dicomPACS®** Upper Cervical Chiropractic Tools. What also impressed him was that **dicomPACS®** could be updated to a full analysis tool for ALL diagnostic imaging needs even MRI or CT images...).

Further information is available under [www.or-technology.com](http://www.or-technology.com)

Doctor Scholten confirmed us numerous advantages of our Upper Cervical software, which we would like to share:

- The ability to follow NUCCA standards of practice in radiographic marking seamlessly and adapt selected points with ease.
- The ability to view and integrate CT and MRI images.
- Direct access to the main server with multiple analysis stations.
- All the extra tools that come with OR Technology's **dicomPACS**<sup>®</sup> software.
- The ability to easily see the raw data and create comparisons.
- The ability to compare multiple images side by side.
- The manner in which the vector calculator works allows easy adaption by the clinician.
- The different and customizable color of the lines make identification clear and easy during a patient's report of findings.
- The speed of service to update the product when issues have occurred has been excellent.

Summarized Jeff Scholten means: "The software's system is modular and so much more than just "another analysis tool". It supports my workflow and makes it easier."

Contact Jeffrey Scholten

[www.vitalposture.com/contact.html](http://www.vitalposture.com/contact.html)

Northland Professional Building

Suite 104, 4600 Crowchild Trail NW

Calgary, T3A 2L6


Phone: 403-247-4257 Fax: 403-247-4275

Email: [info@vitalposture.com](mailto:info@vitalposture.com)

Dr. Scholten on YouTube:

[www.youtube.com/user/calgaryuppercervical](http://www.youtube.com/user/calgaryuppercervical)

Further information is available under  
[www.or-technology.com](http://www.or-technology.com)





Upper Cervical (NUCCA)  
and Standard Chiropractic Tools -  
Diagnostic tools for  
optimal treatment



*dicom* PACS®

The **dicomPACS**® Upper Cervical (NUCCA) Chiropractic tool set has been created in cooperation with leading NUCCA experts from the US and Canada. It offers a variety of ways to reach a fast and accurate diagnosis. Templates like the Cephalometer, Grid, Circumscale, and Relatoscope enable you to continue working as you are used to.

The Chiro Tools offer great possibilities for diagnosing accurately as well as for planning further treatment. According to the tool used, automated center lines and points, defined curves, angle measurements etc., are generated after the manual selection of the points of interest.

 OR Technology

[www.or-technology.com](http://www.or-technology.com)

OR Technology (Oehm und Rehbein GmbH), 18057 Rostock, Germany, Neptunallee 7c  
Tel. +49 381 36 600 500, Fax +49 381 36 600 555  
[www.or-technology.com](http://www.or-technology.com), [info@or-technology.com](mailto:info@or-technology.com)

Info-Hotline: +49 381 36 600 600

